

Ins^{Eye}ts

FAQ

October 2024



What is Cellebrite Inseyets?

Cellebrite Inseyets is the definitive digital forensics suite that includes the industry leading capabilities needed to close more cases.

What does Cellebrite Inseyets include?

Cellebrite Inseyets.UFED

Access devices and data previously unreachable (when compared to UFED) and extract the Full File System data including containerized applications and encrypted data, when compared to UFED.

Cellebrite Inseyets.PA

Process and analyze evidence quickly from the broadest range of apps and sources, enrich the case with in-depth information and validate your evidence easily.

Cellebrite Streamline

Easily and quickly process a device from extraction to UFDR in a few clicks.

Cellebrite Autonomy

Scale-up investigations and run multiple automated processes quickly and accurately at once.
(Please consult with your dedicated sales manager to determine if you qualify for the Autonomy solution as part of your Cellebrite Inseyets purchase)

Cellebrite Quick View

Get a glimpse into the device before extraction and help identify the device's relevance and priority.

Cellebrite Triage

Leverage customizable search criteria to quickly flag relevant devices and significantly improve device prioritization.

Cellebrite PA Cloud

Extract, preserve and analyze cloud-based data using device-based tokens or known credentials and gain access to valuable evidence that is not available on the device.

Cellebrite Commander

Manage permissions and capabilities for all users, including software updates and license updates.

Cellebrite Reader

Review evidence thoroughly, use advanced search, filtering, timeline, entity bookmarks and project saving features to find key details and capture review progress.

How will Cellebrite Inseyets benefit me compared to what I have now?

When upgrading to Cellebrite Inseyets, you will enjoy these added capabilities:

- **UFED customer** – Advanced device access, including unlock-ready, Full File System extractions, Streamlined workflows, Quick Views before full extraction, Triage prioritization, Automated multi-processing capabilities, Cloud-account access, Commander and more.
- **Premium customer** – Streamlined workflows, Quick Views before full extraction, Triage prioritization, Automated multi-processing capabilities, Cloud-account access, Commander and more.

How is this going to affect my current Cellebrite environment with regards to connectivity, compliance, and security requirements?

To take full advantage of Cellebrite Inseyets, a connection to a network is required (there is no exchange of device data). Cellebrite Inseyets is also available as an offline configuration to match your agency's needs. We advise following Cellebrite Inseyets security requirements to prepare your workstation.

What do I need to know when upgrading my existing products to Cellebrite Inseyets?

- **UFED 4PC** – you will require an upgrade kit, and make sure you have Cellebrite Max dongle.
- **Touch 2** - Cellebrite Inseyets isn't supported on Touch 2, we recommend moving to 4PC or computer-based version hardware where you will require a full Cellebrite Inseyets kit, and make sure you have Cellebrite Max dongle.
- **Touch 3** - You will require an upgrade kit, and make sure you have Cellebrite Max dongle (please note that Touch 3 isn't optimal for Cellebrite Inseyets and some capabilities won't work as well).

Is Unlocking included in Inseyets?

Inseyets is unlock-ready, however you need to purchase the unlocks package separately.

What is the Turbo Link?

Turbo Link serves as the advanced connection between the device and the capabilities of the solution. It embodies the cutting-edge technology powering Cellebrite Inseyets, enhancing the speed and efficiency of data extraction and analysis (formerly known as the "Cheetah" adapter).



Does Celebrite Insecrets require new hardware?

Some systems may require additional hardware. Please refer to the product technical sheets for the specific details.

What dongle do I need to be able and upgrade to Celebrite Insecrets?

Celebrite Insecrets requires a Celebrite Max dongle. Check your dongle type [here](#).

What is the difference between PA 7/PA 10 (Celebrite Insecrets.PA)?

Celebrite Insecrets.PA (version 10.x), operates on a database that handles large volumes of data. It introduces new features like Geo-Location Semantics, Media Origin, Import/Export and quicker case reopening.

With the new Celebrite Insecrets offering, you will also enjoy additional enhancements like Cloud and Crypto features. (Celebrite Insecrets.PA and PA 7, share the same decoding abilities)

Can I try Celebrite Insecrets with my current license?

You may explore Celebrite Insecrets.PA v10.x, however, please be aware that you will only be able to utilize the decoding and analysis features. The rest of Celebrite Insecrets capabilities such as extraction and other advanced capabilities will not be available to try.

About Cellebrite

Cellebrite's (Nasdaq: CLBT) mission is to enable its customers to protect and save lives, accelerate justice, and preserve privacy in communities around the world. We are a global leader in Digital Intelligence solutions for the public and private sectors, empowering organizations in mastering the complexities of legally sanctioned digital investigations by streamlining intelligence processes. Trusted by thousands of leading agencies and companies worldwide, Cellebrite's Digital Intelligence platform and solutions transform how customers collect, review, analyze and manage data in legally sanctioned investigations.

- To learn more visit us at www.cellebrite.com
- Contact Cellebrite globally at www.cellebrite.com/contact

